

DORO®

COBRA




BLACK
FOREST
MEDICAL
GROUP

Retractor System
Item no. 1201.100



THE SLIM ALTERNATIVE

CONNECTIVITY GUARANTEED

DORO® Blades
fit all DORO® and
DORO COBRA®
Retractor Arms.

360° BODY ROTATION

offers additional
freedom of
movement, so
the retractor
arm never
interferes with
your procedure.

ONE SIDE, TWO RETRACTORS

The interface of
the DORO COBRA®
Quick-Clamp allows
you to attach up to
two retractor arms
at the same time



REPOSITION WITH EASE

The DORO COBRA®
Flexible Rotary Arms
can be repositioned
anytime without
having to release
the applied
tension.

SNAP-ON ATTACHMENT

With an open screw,
the retractor arm
attaches effortlessly
to find the perfect
position on the
interface

MAXIMISED STABILITY

The formfitting
components of the
system provide a
rigid connection to
the skull clamp, even
over multiple layers
of sterile drape

MIX & MATCH

The available
Mix & Match accessories
allow you to attach
your DORO COBRA®
Retractor Arms rigidly
to the radiolucent
Quick-Rail®, as on the
DORO® Radiolucent
Skull Clamp (3034-00).



DORO COBRA®
Quick-Clamp for
Radiolucent Quick-Rail®
Item no. 1201.022



The retractor system
attaches perfectly to
the Quick-Rail® of the
new DORO LUCENT®
Skull Clamp without any
additional accessories.

DORO COBRA® Retractor System for Aluminum Quick-Rail®

Item no. 1201.100

Included Items	Item no.
DORO COBRA® Quick-Clamp for Aluminum Quick-Rail® (2 pcs)	1201.021
DORO COBRA® Flexible Rotary Arm 420 mm (2 pcs)	1201.008
DORO® Blade, 102 x 7 mm (2 pcs)	3116-111
DORO® Blade, 102 x 10 mm (2 pcs)	3116-121
DORO® Retractor Sterilization Tray	1201.091



BLACK
FOREST
MEDICAL
GROUP

Black Forest Medical GmbH

Bötzingen Straße 86, 79111 Freiburg im Breisgau, Germany
Phone: +49 761 384 222 - 10 | info@blackforestmedical.com

USA Black Forest Medical North America, Inc.
4529 S.E. 16th Place, Suite 101, Cape Coral, FL 33904
Phone: +1 239 369 - 2310 | info.us@blackforestmedical.com

© Black Forest Medical GmbH, 2026 / 01
All rights reserved. The information provided was compiled
with utmost care. Errors can, however, not be excluded.
Errors do not entitle to any legal claims. We reserve
the right to modify the shown products at any time.



Take a peek on our work and visit us on social media!

blackforestmedical.com

Retractor System

Item no. 1201.100

