

DORO COBRA® RETRACTOR SYSTEM



DORO COBRA² Retractor System

Item no. 1201.100



Please find further information on the product's web page.

DORO COBRA[®]
Flexible Rotary Arm
provides full 360° body
rotation while mounted.

Maximized stability
through formfitting
components.
Provides a rigid
connection even over
multiple layers of
sterile drape.

DORO COBRA[®] Flexible Rotary Arms
can be repositioned under tension
without losing the tension.

Snap-on attachment
of **DORO COBRA[®]**
Retractor Arms.

The **DORO COBRA[®]** Quick-Clamp
for Aluminum Quick-Rail[®]
allows the direct attachment of two
DORO COBRA[®] Flexible Rotary
Arms to the Quick-Rail[®] of **DORO[®]**
Aluminum Skull Clamps.



ENHANCE THE VARIABILITY

of your system with the optionally available **DORO COBRA®** Quick-Clamp for Radiolucent Quick-Rail®. This Quick-Clamp alteration allows you to directly connect **DORO COBRA®** Retractor Arms to the Quick-Rail® of the **DORO®** Radiolucent Skull Clamp.

UNFOLD THE FULL POTENTIAL

by combining different elements of the **DORO COBRA®** and the **DORO LUNA®** system.



DORO LUNA® Quick-Clamp
for Radiolucent Quick-Rail®
Item no. 1203.009

DORO COBRA® Quick-Clamp
for Radiolucent Quick-Rail®
Item no. 1201.022



Black Forest Medical GmbH
Headquarters | Bötzingen Straße 86 | 79111 Freiburg, Germany
Phone +49 761 384 222 10 | Email info@blackforestmedical.com

US Subsidiary | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904
Toll Free 877 225 4086 | Email info.us@blackforestmedical.com

Version 03 © Black Forest Medical GmbH, 07/2025,
Freiburg, Germany.

All rights reserved. The information provided was compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims. We reserve the right to modify the products at any time.



blackforestmedical.com