



BLACK
FOREST
MEDICAL
GROUP

DORO COBRA[®] RETRACTOR SYSTEM



DORO COBRA²® Retractor System

Item no. 1201.100



Please find further information on the product's web page.

DORO COBRA[®]

Flexible Rotary Arm provides full 360° body rotation while mounted.



Maximized stability through formfitting components. Provides a rigid connection even over multiple layers of sterile drape.

DORO COBRA® Flexible Rotary Arms can be repositioned under tension without losing the tension.

Snap-on attachment of DORO COBRA® Retractor Arms.

The DORO COBRA® Quick-Clamp for Aluminum Quick-Rail® allows the direct attachment of two DORO COBRA® Flexible Rotary Arms to the Quick-Rail® of DORO® Aluminum Skull Clamps.





ENHANCE THE VARIABILITY

of your system with the optionally available **DORO COBRA®** Quick-Clamp for Radiolucent Quick-Rail®. This Quick-Clamp alteration allows you to directly connect **DORO COBRA®** Retractor Arms to the Quick-Rail® of the **DORO®** Radiolucent Skull Clamp.

UNFOLD THE FULL POTENTIAL

by combining different elements of the **DORO COBRA®** and the **DORO LUNA®** system.



DORO LUNA® Quick-Clamp
for Radiolucent Quick-Rail®
Item no. 1203.009

DORO COBRA® Quick-Clamp
for Radiolucent Quick-Rail®
Item no. 1201.022



pro med instruments – part of Black Forest Medical Group
Headquarters | Bötzingen Straße 86 | 79111 Freiburg, Germany
Phone +49 761 384 222 10 | Email info@blackforestmedical.com

US Subsidiary | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904
Toll Free 877 225 4086 | Email info.us@blackforestmedical.com

Version 03 © 2022 pro med instruments GmbH, Freiburg, Germany.

All rights reserved. The information provided was compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims. We reserve the right to modify the products at any time.



blackforestmedical.com