

DORO®

SKULL PINS



BLACK
FOREST
MEDICAL
GROUP

Find more information at
blackforestmedical.com



DORO® SKULL PINS

DORO® Skull Pins are specifically designed for their respective intended use. Whether it is a conventional surgery or an intra-operative application in conjunction with CT, MR imaging and Fluoroscopy, a matching DORO® Skull Pin is always available.

Within the imaging segment, DORO LUCENT® Skull Pins take the lead, providing a highly visible enhancement compared to the excellent imaging results of the DORO® Titanium Skull Pins.

PRECISELY MACHINED


Easy tip positioning during the pinning procedure

UNIQUE WINGS

Ensure secure fixation in the pin holder

GRIP DESIGN

allows easy placement & removal.



**DORO® Disposable Skull Pins
Stainless Steel Adult**
Color: Blue
Item no. 3006-00
(36 pins), sterile



**DORO® Disposable Skull Pins
Stainless Steel Pediatric**
Color: Yellow
Item no. 3006-10
(36 pins), sterile



**DORO® Disposable Skull Pins
Stainless Steel Adult**
Color: Black
Item no. 3006-50
(36 pins), sterile



**DORO® Reusable Skull Pins
Stainless Steel (3 Pins)**
Ring-Color: Black
Item no. 3005-00 (Adult)
Item no. 3004-00 (Pediatric)

NON-IMAGING



**DORO® Disposable Skull Pins
Titanium Adult**
Color: Orange
Item no. 3006-20
(36 pins), sterile



**DORO® Disposable Skull Pins
Titanium Pediatric**
Color: Green
Item no. 3006-30
(36 pins), sterile



**DORO LUCENT® Disposable Skull Pins
Adult**
Color: Black
Item no. 1106.003 (3 pins)*
Item no. 1106.103 (9 pins)



**DORO® Reusable Skull Pins
Titanium (3 Pins)**
Ring-Color: Blue
Item no. 3005-50 (Adult)
Item no. 3004-50 (Pediatric)

IMAGING

*Not available in the United States & Canada



Black Forest Medical GmbH
Böttinger Straße 86, 79111 Freiburg, Germany
Phone: +49 761 384 222 - 10 | info@blackforestmedical.com

USA Black Forest Medical North America, Inc.
4529 S.E. 16th Place, Suite 101, Cape Coral, FL 33904
Phone: +1 239 369 - 2310 | info.us@blackforestmedical.com

© Black Forest Medical GmbH, 2025 / 09
All rights reserved. The information provided was compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims. We reserve the right to modify the shown products at any time.



blackforestmedical.com