

DORO[®]

LUCENT



Cranial Stabilization Set IMRIS

Item no. 4003.600



BLACK
FOREST
MEDICAL
GROUP

EFFICIENT WORKFLOWS

The DORO LUCENT® iMRI Cranial Stabilization Set IMRIS is a lightweight and intuitive intraoperative solution for the IMRIS Invision™ Surgical Theatre, enabling clear and accurate visualizations for precise diagnoses and treatment. The interface to the IMRIS coil can be adjusted continuously, allowing for precise positioning. The system is fast and easy to set-up and its DORO® exclusive Quick-Rail® facilitates the easy adaptation of accessories.

The DORO LUCENT® Headplate allows the clinical staff to position the patient as required for surgery before starting the actual pinning procedure. This provides an optimized OR workflow and less stress for the patient.

Automatic Image Registration Ready

The system is compatible* with image-guided surgery (IGS) systems from the leading companies in the field, including the Automatic Image Registration (AIR)** technology from Brainlab. It allows neurosurgeons to seamlessly integrate reacquired intra-operative MR images, eliminating the need to access anatomical landmarks.

* Compatibility subject to the type of IMRIS MRI scanner used.
** Automatic Registration iMRI and accessories.

INCLUDED ITEMS	
DORO® Skull Clamp Radiolucent for iMRI	4003.001
DORO LUCENT® Headplate	4003.005
DORO® Adjusting Wrench	4003.011
DORO LUCENT® Storage Case	4003.030
DORO LUCENT® iMRI Interface Adaptor IMRIS	4003.611
DORO LUCENT® iMRI Coil Connector IMRIS	4003.631

COMPATIBLE* WITH AUTOMATIC IMAGE REGISTRATION**

Fast integration of MR images without the need to access anatomical landmarks



PLAY-FREE OPEN/LOCK

Enhanced precision for neuro-navigation

IMPROVED FLEXIBILITY

Slim and thoughtful design allows for easy and collision-free accessibility

UNIQUE QUICK-RAIL®

Intuitive mounting of retractor systems or other accessories



HEADPLATE INCLUDED

Optimized OR workflow and reduced stress for the patient

OPTIMIZED FOR IMRI

Continuous adjustability for precise adjustment of coil position

Minimized amount of titanium results in less artifacts.

DORO® LUCENT



A NEW ERA OF IMAGING QUALITY

The combination of the DORO LUCENT® iMRI Cranial Stabilization Set IMRIS and its accompanying DORO LUCENT® Skull Pins unveils an unparalleled level of detail and accuracy, providing exceptional support to neurosurgeons and their teams during intra-operative MR imaging procedures. When compared to other products on the market, the remarkable clarity and near elimination of imaging artifact truly sets it apart.





BLACK
FOREST
MEDICAL
GROUP

Black Forest Medical GmbH

Bötzingen Straße 86, 79111 Freiburg, Germany
Phone: +49 761 384 222 - 10 | info@blackforestmedical.com

USA Black Forest Medical North America, Inc.
4529 S.E. 16th Place, Suite 101, Cape Coral, FL 33904
Phone: +1 239 369 - 2310 | info.us@blackforestmedical.com

© Black Forest Medical GmbH, 2025 / 09
All rights reserved. The information provided was compiled
with utmost care. Errors can, however, not be excluded.
Errors do not entitle to any legal claims. We reserve
the right to modify the shown products at any time.



Take a peek on our work and visit us on social media!

blackforestmedical.com

Cranial Stabilization Set IMRIS

Item no. 4003.600

