

# DORO® NAVIGATION ADAPTORS

To ensure a secure and accurate navigation-assisted operation, navigation systems are dependent on important requirements:

- They have to be stable, so that a precise reference point is guaranteed but flexibly adjustable as well.
- The reference frame must be high enough to be captured by the camera and simultaneously close to the surgical site.
- Also they must be small-sized to allow intra-op imaging.

**DORO® Easy Connect Navigation Adaptors** and **DORO® Navigation Adaptors Radiolucent** are designed to comply with all these requirements in a very effective manner.

## DORO® NAVIGATION ADAPTORS

- Connect directly to the integrated navigation starburst on the Quick-Rail® of the **DORO® QR3 Skull Clamp**, **DORO® Multi-Purpose Skull Clamp** and **DORO® Skull Clamp Teflon®**.
- Provide a very stable and secure connection of the navigation interface.
- Allow to position the navigation interface close to the surgical target.





Please find further information on the product's web page.

# DORO® NAVIGATION READY

## DORO® EASY CONNECT NAVIGATION ADAPTORS

Adaptors for the  
DORO® Aluminum, Multi-Purpose and Teflon® Skull Clamps.

For navigation where intra-operative imaging is not required – the 1204 series.



DORO® Easy Connect  
Navigation Adaptor Stryker®  
Item no. 1204.001



DORO® Easy Connect  
Navigation Adaptor Brainlab  
Item no. 1204.002



DORO® Easy Connect  
Navigation Adaptor Medtronic  
Item no. 1204.003

## DORO LUCENT® NAVIGATION ADAPTORS

Radiolucent adaptors for the DORO® Radiolucent Skull Clamp for intra-operative imaging applications – the 3033 series.



DORO LUCENT®  
Navigation Adaptor, Stryker®  
Item no. 3033-64



DORO LUCENT®  
Navigation Adaptor, Medtronic  
Item no. 3033-65



DORO LUCENT®  
Navigation Adaptor, Brainlab  
Item no. 3033-66  
for Brainlab Reference Array  
19145-01 with or without  
elongation 19145-02



pro med instruments – part of Black Forest Medical Group  
Headquarters | Bötzingen Straße 86 | 79111 Freiburg, Germany  
Phone +49 761 384 222 10 | Email info@blackforestmedical.com

US Subsidiary | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904  
Toll Free 877 225 4086 | Email info.us@blackforestmedical.com

Version 04 © 2021 pro med instruments GmbH, Freiburg, Germany.

All rights reserved. The information provided was compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims. We reserve the right to modify the products at any time. The present flyer describes products which may not be approved for sale in all countries.



blackforestmedical.com