



# DORO® QR3 SKULL CLAMP



# DORO® QRR3 Skull Clamp



**EASY TO READ**  
laser etched scale  
(N and lbs) on the  
torque screw.




**NAVIGATION READY**  
Stainless steel starbursts  
on both sides of the Quick-  
Rail® allow direct attachment  
of movement-sensitive  
navigation tracking device  
adaptors.



**QUICK-RAIL®**  
An easy to locate, easy to handle  
interface for a stable and secure  
connection of **DORO LUNA®**  
and **DORO COBRA®** Retractor  
Systems.





**ULTRA-PRECISE  
FABRICATION**

No play of two pin rocker arm when locked.

**LARGE DORO® SKULL PIN  
PORTFOLIO**

Disposable and reusable skull pins for adult and pediatric applications are available.

**ULTRA-PRECISE  
FABRICATION**

Minimum play between skull clamp and extension arm.

**DORO® QR3 Skull Clamp**  
item no. 1001.001



**MIX & MATCH**

metal and radiolucent parts by means of a **DORO® Interface Adaptor**, item no. 3033-51.



Please find further information on the product's web page.



The **DORO®** QR3 Skull Clamp provides 3-pin rigid cranial fixation and features the latest generation of Quick-Rail®.

Designed with straight and upright proportions, the QR3 offers more Quick-Rail® space than all its predecessors.

Furthermore, the improved Quick-Rail® design makes it easier to connect accessories, even through multiple layers of sterile drape.



pro med instruments – part of Black Forest Medical Group  
Headquarters | Bötzingen Straße 86 | 79111 Freiburg, Germany  
Phone +49 761 384 222 10 | Email [info@blackforestmedical.com](mailto:info@blackforestmedical.com)

US Subsidiary | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904  
Toll Free 877 225 4086 | Email [info.us@blackforestmedical.com](mailto:info.us@blackforestmedical.com)

Version 04 © 2021 pro med instruments GmbH, Freiburg, Germany.

All rights reserved. The information provided was compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims. We reserve the right to modify the products at any time.



[blackforestmedical.com](http://blackforestmedical.com)