



BLACK
FOREST
MEDICAL
GROUP

DORO®

QUICK-CLAMP WITH STARBURST

- FOR RADIOLOGENT QUICK-RAII® -



DORO® QUICK-CLAMP WITH STARBURST

- FOR RADIOCULENT QUICK-RAIL® -

The DORO® Quick-Clamp with Starburst for Radiolucent Quick-Rail® serves as an intuitive interface* for devices used in frameless stereotactic or endoscopic procedures.

Designed for use over the sterile barrier and formfittingly connected to the proven Quick-Rail® of the DORO® Skull Clamp Radiolucent, the adaptor provides a stable and secure connection point for third party devices.



Please find further information on the product's web page.



DIRECTLY CONNECT stereotactic devices to the Quick-Rail®

*The interface is designed to bear a max. load of 16.8 Nm (which equals a weight of 3.4 kg attached to a horizontally mounted, 50 cm long arm). It is validated for use with the Brainlab VarioGuide™. It is the responsibility of the user to verify that third party devices can safely be secured to the DORO® Quick-Clamp with Starburst for Radiolucent Quick-Rail®. Product is shown with Brainlab VarioGuide™ and DORO® Skull Clamp Radiolucent.



pro med instruments – part of Black Forest Medical Group
Headquarters | Bötzingen Straße 86 | 79111 Freiburg, Germany
Phone +49 761 384 222 10 | Email info@blackforestmedical.com

US Subsidiary | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904
Toll Free 877 225 4086 | Email info.us@blackforestmedical.com

Version 04 © 2022 pro med instruments GmbH, Freiburg, Germany.

All rights reserved. The information provided was compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims. We reserve the right to modify the products at any time. VarioGuide™ is a trademark of Brainlab®.