

DORO®



QR3 Headrest System

Item no. 1001.100



BLACK
FOREST
MEDICAL
GROUP



PRECISE
NAVIGATION

Starbursts on both sides of the Quick-Rail® allow direct attachment of navigation tracking device adaptors.



OPTIONAL
INTERFACE

Utilizing the DORO® Side Rail Adaptor, the system can be effortlessly mounted to OR table side rails.

EFFORTLESS
CLEANING

Dedicated drain holes on the base unit allow for better flushing and drainage.



PLAY-FREE
OPEN/LOCK

Enhanced precision for neuro-navigation.

MORE SPACE

Unique shape ensures sufficient free space to place the intubation tube, even in prone position.

REFINED
APPEARANCE

Slim and thoughtful design enables intuitive handling and provides collision-free accessibility

EXCELLENT
DURABILITY

High quality material and process standards result in increased resistance to operational stress.

ALWAYS
RELIABLE

Designed for daily use, the DORO® QR3 Headrest System streamlines all non-imaging neurosurgical procedures. Fast and easy to set up, this system will be your everyday companion. Deliberately designed features like built-in navigation adaptor interfaces or its Mix & Match compatibility with DORO® Radiolucent products make it an incredibly versatile device. Safe for your patients, convenient for you.

A Proven Skull Clamp

The DORO® QR3 Skull Clamp provides 3-pin rigid cranial fixation and features the latest generation of Quick-Rail®. Designed with straight and upright proportions, the QR3 offers more Quick-Rail® space than all its predecessors.

INCLUDED ITEMS	
DORO® QR3 Skull Clamp	1001.001
DORO® Swivel Adaptor Aluminum	3002-00
DORO® Base Unit Aluminum	3001-00

OPTIONALLY AVAILABLE	
DORO® Disposable Skull Pins Stainless Steel, Adult	3006-00
DORO® Disposable Skull Pins Stainless Steel, Pediatric	3006-10
DORO® Swivel Horseshoe Headrest System, Adult	3009-00
DORO® Swivel Horseshoe Headrest System, Pediatric	3008-00
DORO® Side Rail OR Table Adaptor	3011-00
DORO® Side Rail OR Table Adaptor AMSCO®	3011-10
DORO® Storage Case Headrest System	1001.030
DORO® Easy Connect Navigation Adaptor Brainlab	1204.001
DORO® Easy Connect Navigation Adaptor Stryker®	1204.002
DORO® Easy Connect Navigation Adaptor Medtronic	1204.003



QUICK-RAIL®

AN INNOVATIVE INTERFACE

The term QR3 refers to the 3rd generation of our innovative Quick-Rail® technology, which is directly implemented in both extension arm and fixed arm of our aluminum skull clamps. Due to its unique shape, it can easily be located, even when covered by multiple layers of sterile drape. The Quick-Rail® serves as a form-fitting, intuitive interface, ensuring reliable attachment of DORO® Retractor Systems and other accessories.



OPTIONAL
HORSESHOE


The non-invasive solution for adult and pediatric treatment.



BLACK
FOREST
MEDICAL
GROUP

Black Forest Medical GmbH

Bötzingen Straße 86, 79111 Freiburg im Breisgau, Germany
Phone: +49 761 384 222 - 10 | info@blackforestmedical.com

 **USA** Black Forest Medical North America, Inc.
4529 S.E. 16th Place, Suite 101, Cape Coral, FL 33904
Phone: +1 239 369 - 2310 | info.us@blackforestmedical.com

© Black Forest Medical GmbH, 09/2025
All rights reserved. The information provided was compiled
with utmost care. Errors can, however, not be excluded.
Errors do not entitle to any legal claims. We reserve
the right to modify the shown products at any time.



Take a peek on our work and visit us on social media!

blackforestmedical.com

QR3 Headrest System

Item no. 1001.100

