

DORO® QR3 HEADREST SYSTEM TEFLON®



DORO® QR3 HEADREST SYSTEM TEFLON®

The DORO® Teflon® Headrest System is a mechanical support system which is used in head and neck surgery. This system allows the patient's head to be positioned and secured during the operation. It can be used for cranial fixation of the patient in diverse positions.

By adapting the Teflon® Headrest System to the QR3 format, all the beneficial features of this state-of-the-art skull clamp design become available for the Teflon® series.

Furthermore the PMI exclusive Teflon® coating allows it to be reprocessed in a washer/disinfector. The complete headrest system conforms to WHO Infection Control Guidelines, is machine washable, disinfectable and autoclavable.

DORO® QR3
HEADREST
SYSTEM TEFLON®
IS AUTOCLAVABLE.




Please find further information on the product's web page.

FEATURES

- New QR3 Quick-Rail® technology for quick connection of accessories like the new DORO LUNA® Retractor System as well as Retractor Arms, J-Arms and Halo Retractor Systems.
 - Autoclavable.
 - Machine washable.
 - Conforms to WHO Infection Control Guidelines for cleaning and sterilization.
 - Navigation ready.
 - 3-pin rigid cranial fixation.
 - For supine, prone, lateral or sitting positions.
- For lateral positioning, we recommend the DORO® Base Unit Parkbench (not Teflon® coated) (item no. 3001-006).
- Sitting positions: Only in conjunction with the DORO® Crossbar Adaptor (not Teflon® coated) (item no. 3007-00).
- Excellent stability and flexibility.
 - Exclusive – only available from PMI.



SYSTEM COMPONENTS



DORO® QR3
Headrest System Teflon®
Item no. 3003-100

DORO® Skull Clamp Teflon®
Item no. 3003-009

DORO® Swivel Adaptor Teflon®
Item no. 3002-009

DORO® Adjustable Base Unit Teflon®
Item no. 3001-009



DORO® Storage Case Headrest System
(optionally available)
Item no. 1001.030



Black Forest Medical GmbH
Headquarters | Bötzingen Straße 86 | 79111 Freiburg, Germany
Phone +49 761 384 222 10 | Email info@blackforestmedical.com

US Subsidiary | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904
Toll Free 877 225 4086 | Email info.us@blackforestmedical.com

Version 03 © Black Forest Medical GmbH, 07/2025, Freiburg, Germany.

All rights reserved. The information provided was compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims. We reserve the right to modify the products at any