



BLACK
FOREST
MEDICAL
GROUP

DORO®

SKULL CLAMP FOR SPINAL FRAMES



DORO® SKULL CLAMP FOR SPINAL FRAMES

– RADIOLUCENT –

This modification of our proven, metal free DORO® Skull Clamp Radiolucent is designed for spinal surgeries and can be mounted to all established spinal OR tables and spinal frames. The new, slim design allows the skull clamp to be lowered between the rails of the open frame, simplifying patient positioning and imaging procedures.

DORO® Skull Clamp Radiolucent
for Spinal Frames
Item no. 4004.101



Please find further information on the product's web page.

MORE FLEXIBILITY IN
PATIENT POSITIONING



Standard DORO® Skull Clamp
Radiolucent mounted to a MIZUHO
OSI Insite table.



DORO® Skull Clamp Radiolucent for
Spinal Frames mounted to a MIZUHO
OSI Insite table.



pro med instruments – part of Black Forest Medical Group
Headquarters | Bötzingen Straße 86 | 79111 Freiburg, Germany
Phone +49 761 384 222 10 | Email info@blackforestmedical.com

US Subsidiary | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904
Toll Free 877 225 4086 | Email info.us@blackforestmedical.com

Version 03 © 2021 pro med instruments GmbH, Freiburg, Germany.
All rights reserved. The information provided was compiled with utmost care. Errors can, however, not be excluded. Errors do not entitle to any legal claims. We reserve the right to modify the products at any time.