

To whom it may concern

### DORO® SKULL PIN compatibility with DORO® products and third-party products

DORO® PRODUCTS are manufactured with high precision and under strict quality requirements. In order to ensure optimum functionality of our system, the various components perfectly match with each other. If the components do not match with each other, safe utilization of the products is not ensured.

DORO® SKULL PINS in particular are manufactured with highest precision to perfectly match the dimensions of our pin receptacles. This ensures easy and safe mounting as well as a secure fit of the DORO® SKULL PINS and prevents them from accidentally falling out. As Black Forest Medical Group, namely pro med instruments GmbH part of Black Forest Medical Group (hereinafter BFMG) as the manufacturer, checks the compatibility of its components continuously, BFMG can ensure the compatibility of all DORO® SKULL PINS with all DORO® SKULL CLAMPS.

The pin receptacle of skull clamps from other manufacturers are usually designed similarly, which means that DORO® SKULL PINS might also be used safely with these products if compatibility is ensured. Such compatibility can only be ensured by BFMG if an individual check of the accurate fit of the DORO® SKULL PINS in pin receptacles of the third-party product has been carried out with a calibrated limit gauge provided by BFMG and used in accordance with these instructions.

Only if the qualified limit gauge confirms the accurate fit of the DORO® SKULL PINS, BFMG can declare product compatibility with regard to DORO® SKULL PINS in connection with the third-party product.

Sincerely,



Matthias Schüle  
Owner/Managing Director



Irina Bauer  
Director Sales International & Marketing

pro med instruments GmbH  
part of Black Forest Medical Group