




# Which DORO LUCENT® iMRI Headrest System works for you?

Check the type of MRI scanner and OR table in your OR environment and compare it to the overview below to identify the DORO LUCENT® iMRI Headrest System that will fit your OR environment.

MRI SCANNER	PHILIPS Ingenia MR-OR (1.5T & 3T)	GE 3.0T Signa HDxt / Signa Artist & Architect / Optima MR450w / Discovery MR750 & 750w	SIEMENS MAGNETOM Aera / Skyra Sola / Vida	SIEMENS MAGNETOM Aera / Skyra	
OR TABLE with Transfer Table Top	MIZUHO	MAQUET / GETINGE Alphaquet & Otesus / Magnus	TRUMPF TruSystem		
					
<b>DORO LUCENT®</b> iMRI Headrest System MIZUHO Item no. 4003.350	<b>DORO LUCENT®</b> iMRI Headrest System MAQUET Item no. 4003.300	<b>DORO LUCENT®</b> iMRI Headrest System TRUMPF Item no. 4003.200			
Included Items	Item no.	Included Items	Item no.	Included Items	Item no.
DORO® Skull Clamp Radiolucent for DORO LUCENT® Parallelogram Adaptor	4003.001	DORO® Skull Clamp Radiolucent for DORO LUCENT® Parallelogram Adaptor	4003.001	DORO® Skull Clamp Radiolucent for DORO LUCENT® Parallelogram Adaptor	4003.001
DORO® Headplate	4003.005	DORO® Headplate	4003.005	DORO® Headplate	4003.005
DORO® Adjusting Wrench for DORO LUCENT® Parallelogram Adaptor	4003.011	DORO® Adjusting Wrench for DORO LUCENT® Parallelogram Adaptor	4003.011	DORO® Adjusting Wrench for DORO LUCENT® Parallelogram Adaptor	4003.011
DORO LUCENT® Storage Case for Parallelogram Adaptors	4003.030	DORO LUCENT® Storage Case for Parallelogram Adaptors	4003.030	DORO LUCENT® Storage Case for Parallelogram Adaptors	4003.030
DORO LUCENT® iMRI Parallelogram Adaptor MAQUET	4003.311	DORO LUCENT® iMRI Parallelogram Adaptor MAQUET	4003.311	DORO LUCENT® iMRI Parallelogram Adaptor TRUMPF	4003.211
DORO LUCENT® Disposable Skull Pins, Adult (3 pcs)	1106.003	DORO LUCENT® iMRI Transfer Collision Indicator	4003.XXX	DORO LUCENT® iMRI Pad TRUMPF	4003.212
		DORO LUCENT® Disposable Skull Pins, Adult (3 pcs)	1106.003	DORO LUCENT® iMRI Transfer Collision Indicator	4003.XXX
				DORO LUCENT® Disposable Skull Pins, Adult (3 pcs)	1106.003
				DORO® Torx Screwdriver T45, nonmagnetic	3033-68

A fitting iMRI Transfer Collision Indicator must be purchased from Phillips.

**PLEASE NOTE:**  
With the first installation of a DORO LUCENT® iMRI Headrest System, an expert training is mandatory (page 32). Please contact your local DORO® representative for details.