


UNLIMITED POSSIBILITIES



DORO LUNA® and DORO COBRA® Retractor Systems can be **directly mounted to the DORO LUCENT® Skull Clamp** (1101.001).

This brings all the advantages of our flagship retractor systems to imaging-supported procedures. When intraoperative imaging is required, the entire system can be temporarily removed quickly and easily by loosening just two thumbscrews.

After imaging is completed, it can be re-mounted in seconds using the same simple process. This allows surgeons to work with a single retractor system, maintaining their familiar routine and workflow from non-imaging procedures.

MIX & MATCH



For the DORO® Skull Clamp Radiolucent (3034-00) and the DORO® Skull Clamp for Spinal Frames (4004.101), additional adaptors are available to correctly mount the retractor system. These adaptors also make it possible to mount your DORO COBRA® Retractor Arms directly to a DORO® Radiolucent horseshoe headrest (e.g. 3035-00):



DORO LUNA®
Quick-Clamp for
Radiolucent Quick-Rail®
Item no. 1203.009



DORO COBRA®
Quick-Clamp for
Radiolucent Quick-Rail®
Item no. 1201.022